

Top Dog Manufacturing

25055 Route 2
Kensington, PE
Canada C0B 1M0



TOLL FREE 1-800-472-2412
EMAIL topdog@topdogmfg.com

TECHNICAL DATA SHEET

Medical-Reusable Level II Isolation Gown - L2IG



Health/Medical
industry

Reusable

Recyclable

Easily
donned/doffed

Material

Care instructions

Colours

Packaging

Size

MADE IN CANADA.

Material and design tested to AAMI PB70 standard. The lightweight, supple material and wrap-round design provide outstanding comfort and protection, ensuring that the wearer remains safe while working in healthcare applications. Cotton tie straps at the back of the neck and waist are generous in width and length for easy tying, donning, and doffing. Choice of thumb loop or elastic cuff wrist band.

For use in health and medical facilities, where high-quality, reusable, and laundry-safe product is needed. Not for sale In Canada

REUSABLE/Launderable. Multiple uses results in lower garment costs and reduced environmental impact.

100% recyclable at end of use cycle. Lower costs; reduced environmental impact.

Donned and doffed by the wearer without touching the possible contact area, as per infection control guidelines. Thumb loop makes for easy donning and doffing of gloves.

100% 70gsm polyester fabric with a 0.02 TPU laminate coating. 10 Mil 250 Microns overall thickness.

This product is a 100% barrier to fluids. 100% unaffected by disinfectants, alcohol-based products, soaps, detergents, and water. Highly resistant to wide range of cleaning chemicals.

Tested to AAMI PB70 standard.

Machine wash in cold or hot water using any disinfectant or soap (CDC recommends oxygen-based detergent at 71C/100F for 25 minutes). Do not use bleach.

Tumble dry to 15-25 minutes, maximum temperature of 80°C/176°F, then 10 minutes low heat tumble.

Yellow.

8 per box; 4 kg (8.8 lbs) per box; bulk pack.

One universal size. 115cm (46") length. Check with us for custom sizing.

